MicroChem Laboratories (Pvt.) Ltd.
No.112/1, 1/1, Stanley Thilakarathne Mawatha
, Nugegoda.

INSTRUCTION FOR SAMPLE COLLECTION FOR MICROBIOLOGICAL ANALYSIS

1. GENERAL

Sample for microbiological analysis must only be collected in sterile bottles or bags. Collect samples that are representative of the water /food being tested, and use aseptic techniques to avoid sample contamination. When sample is collected, leave ample air space in the bottle (at least 2.5 cm), and an identification number shall be marked on each container.

1.1 ORDER OF TAKING SAMPLE

When a number of samples for various purposes are being collected from same sampling point, the sample for **microbiological examination should be collected first.**

2. WATER

2.1 SIZE OF THE SAMPLE

The volume of a sample shall be at least 3/4 of the bottle for water & beverage, 250g for food items.

2.2 OPENING AND FILLING THE BOTTLE

The bottles shall be kept unopened until moment it is required for filing. During sampling, the bottle shall be held by the base in one hand, while the other hand the stopper and cover are removed together, The bottle shall be filled, without rinsing, and the stopper shall be replaced immediately.

2.3 TRANSPORT AND STORAGE OF SAMPLE

Examination of bacteriological quality shall preferably be started within 1 hour of Collection of sample. In any case the interval between collection of the sample and beginning of the examination shall not be exceed 24 hours.

If the time interval between collection and examination is to exceed 1 hour the sample shall be transported to laboratory in iced cooler (ideally sample should be held Less than 6 °C). Sample shall not be frozen.

3. FOOD & OTHER

3.1 SAMPLE COLLECTION

- | | |
|---------------------------------|---|
| Finish products | - Consumer package of foods should be sampled from original previously unopened container. |
| Bulk liquid material | - Aseptically transfer a representative sample portion to a sterile container. |
| Bulk solid or semi solid | - Aliquots from several areas of the food under examination should be taken to ensure a representative samples sterile spoon or spatula should be used. |
| Frozen bulk material | - Obtain sample from few areas of the food with a sterile sharp sampling Instrument. |

3.2 TRANSPORT & STORAGE

Sample should be delivered to the laboratory as soon as possible. Refrigerated products must be transported & stored in refrigerated condition 4 ± 2 °C and frozen product be transported and stored in frozen state.

N.B.

For Chemical analysis , water samples should be provided in a separate bottle. Volume should not be less than 1L.

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